



LiNX Output Module



There is no Output Module for LINX listed within the chair **Setup Menu** like every other module listed.

TO BEGIN:

Make a Function Card that is a Utility (Purple) Card from the menu once the Programming tool is unlocked. Write to the chair.

Once rebooted and ready, go into the **Utility Function Card** you just created and open the card.

THE FOLLOWING INTERFACE WILL BE SHOWN:

- Utility Function Name:** This will be the name of the Function Card just created.
- Enable Utility Function:** This setting will enable the card within the profile.
- Function User Input:** This setting determines what Driver Control will control this function.
- Navigation Timeout Enabled:** Default is ON. When enabled, this sets the card to timeout to the menu.
- Allow Multiple Quadrants:** Default is OFF. When turned ON, this setting will allow quadrants to be activated simultaneously.

The next section contains the programming for the 4 quadrants: **FORWARD, REVERSE, RIGHT AND LEFT.**

** PLEASE NOTE: If you are creating a card for anything other than Lights or Horn, you will need to set the Outputs.*

Forward Momentary / Short Press	
Output:	Output 1 (FORWARD)
Activation Mode:	Momentary
Display Icon:	Choose which icon you want to represent the action

Reverse Momentary / Short Press	
Output:	Output 2 (REVERSE)
Activation Mode:	Momentary
Display Icon:	Choose which icon you want to represent the action

Left Momentary / Short Press	
Output:	Output 3 (LEFT)
Activation Mode:	Momentary
Display Icon:	Choose which icon you want to represent the action

Right Momentary / Short Press	
Output:	Output 4 (RIGHT)
Activation Mode:	Momentary
Display Icon:	Choose which icon you want to represent the action

LINX UTILITY FUNCTION CARD SETTINGS

QUADRANT	OUTPUT NUMBER	ACTIVATION TYPE
● Forward Quadrant	● Output 1	● Momentary
● Reverse Quadrant	● Output 2	● Momentary
● Left Quadrant	● Output 3	● Momentary
● Right Quadrant	● Output 4	● Momentary





ASL Wireless Mouse Emulator Setup for LiNX Module




THREE (3) SWITCH: DIP SWITCHES 11 and 12

* This configuration is typically used with a Head Array.
There is NO Reverse Pad so the Reverse cable is not used.



Cable Jacks / Plugs	Port on Mouse Emulator	Behavior
Forward	Forward	Forward / Reverse Toggle
Reverse	N/A	
Left	Left Click	Left Click / Double Click
Right	Right	Right / Left Toggle
Additional Switch		

Diagram of ASL Wireless Mouse Emulator
Dip Switch Panel



FOUR (4) SWITCH: DIP SWITCHES 11 and 12

* This configuration is used with a joystick whenever a client cannot let go of a joystick for clicks.



Cable Jacks / Plugs	Port on Mouse Emulator	Behavior
Forward	Forward	Forward / Reverse
Reverse	Left Click	Left Click / Double Click
Left	Left	Left
Right	Right	Right
Additional Switch		

FOUR (4) SWITCH: DIP SWITCHES 11 and 12

* This configuration is used with a joystick when a client is able to let go of a joystick and wants to use an external switch for clicking.



Cable Jacks / Plugs	Port on Mouse Emulator	Behavior
Forward	Forward	Forward / Reverse
Reverse	Left Click	Left Click / Double Click
Left	Left	Left
Right	Right	Right
Additional Switch	Right Click	Right Click



ASL Wireless Mouse Emulator Setup for LiNX Module




FIVE (5) SWITCH: DIP SWITCHES 11 and 12



* This configuration is selected for clients using a joystick and would like one or both switches configured for LEFT, RIGHT, and Double Clicking.

Diagram of ASL Wireless Mouse Emulator Dip Switch Panel



Cable Jacks / Plugs	Port on Mouse Emulator	Behavior
Forward	Forward	Forward
Reverse	Reverse	Reverse
Left	Left	Left
Right	Right	Right
	Right Click	Right Click
Additional Switch	Left Click	Left Click

	3 SWITCH	4 SWITCH	4 SWITCH Alternative	5 SWITCH
DIP SWITCHES 11 & 12				
802.4 Cable Forward	Forward Jack (Toggles Up / Down)	Forward Jack (Toggles Up / Down)	Forward Jack (Toggles Up / Down)	Forward Jack
802.4 Cable Reverse		Left Click Jack	Left Click Jack	Reverse Jack
802.4 Cable Left	Left Click Jack	Left Jack (Toggles Up / Dn)	Left Jack	Left Jack
802.4 Cable Right	Right Jack (Toggles Right / Left)	Right Click Jack	Right Jack	Right Jack
Additional Switch			Right Click Jack	Left Click Jack
Additional Switch				Right Click Jack
	<i>Used primarily with Head Arrays. There is NO Reverse Pad, so the Reverse Plug is NOT used.</i>	<i>Used with a Joystick when a client is unable to let go of a joystick for Clicks.</i>	<i>Used with a Joystick when a client is able to let go of joystick and would like to use an External Switch for Clicking.</i>	<i>Suggested for clients using a Joystick and would like One or Both Switches configured to Left, Right, or Double Clicking.</i>