



# ASL121S **ATOM** ALS PACKAGE

EXPANDABLE DRIVE CONTROL SYSTEM



**ASL 121S ATOM**  
ALS Package on Chair



**ASL 121S ATOM**  
ALS Drive Control System



## DIGITAL DRIVE RESET SWITCHES

## PRODUCT SPECIFICATIONS

The ASL 121 ATOM ALS package is a digital head control system specifically designed for the need of a person diagnosed with ALS.

This deconstructed ATOM Head Array allows a person to rest and be supported on the back pad and the forward drive direction can be mounted at the knee or hands depending on the most consistent needs and location for the patient.

This system will allow a person to control their chair fully and functionally. The adjustability options are designed to fit customizable needs and progressions. The ATOM ALS Head Array is designed for persons with ALS or diagnoses that have quickly changing needs.

With the addition of the ASL accessories, a client can wirelessly connect to an ACC device, mouse emulator and the Tecla E. device, mouse emulator and the Tecla E.

<p><b>Pulse Mechanical Switch</b></p> <p>ASL 320</p>	<p><b>Ultra Light Switch</b></p> <p>ASL 314</p>	<p><b>Non-Adjustable Proximity Sensor</b></p> <p>ASL 204</p>	<p><b>Wobble Switch Light Pressure</b></p> <p>ASL 304A</p>
--	---	--	--

## MOUNTING HARDWARE

<p><b>Multi Axis Headrest Mounting Hardware</b></p> <p>ASL 606</p>	<p><b>Arm Rest Switch Mounting Bracket</b></p> <p>ASL 611</p>	<p><b>Ultra Sub-Occipital</b></p> <p>ASL 908</p>	<p><b>Swing Away Switch Mounting Assembly - Left</b></p> <p>ASL 930L</p>
<p><b>Swing Away Switch Mounting Assembly - Right</b></p> <p>ASL 930R</p>	<p><b>Swing Away Switch Mounting Assembly - Both</b></p> <p>ASL 930B</p>		

## ATOM - ACCESSORIES

<p><b>ATOM Wireless Triple Switch Receiver</b></p> <p>ASL 557-3</p>	<p><b>ATOM Wireless Mouse Emulator</b></p> <p>ASL 558</p>	<p><b>Tecla E Bluetooth Interface</b></p> <p>ASL 545E</p>
---	---	---